

Aquagent® Titrant 2, free from pyridine, titrant-component for volumetric Karl Fischer titration**Identification**

Taric code: 3822 00 00

Applications

analytical chemistry, laboratory reagent.

Specifications

For use with:

Aquagent® Solvent (AQ0002)

Contains methanol

Physical data

- Flash point: 11 °C
- Hygroscopic

Safety - GHS

Signal Word: Danger

Hazard Statements:

H225: Highly flammable liquid and vapour.

H331: Toxic if inhaled.

H370: Causes damage to organs.

**Precautionary Statements:**

P210: Keep away from heat / sparks / open flames / hot surfaces. - No smoking.

P241: Use explosion-proof electrical / ventilating / lighting / equipment.

P260: Do not breathe dust / fume / gas / mist / vapours / spray.

P303+P361+P353: IF ON SKIN (or hair): Remove / Take off immediately all contaminated clothing. Rinse skin with water / shower.

P405: Store locked up.

P501a: Dispose of contents / container in accordance with local / regional / national / international regulations.

Toxicological data

- LD 50 (oral, rat): 5628 mg/kg (methanol)

Transport/storage

- ADR: 3 FT1 II • UN 1992 • FLAMMABLE LIQUID, TOXIC, N.O.S.(Aquagent® Titrant 2, titrant-component for Karl Fischer titration)
- IMDG: 3 II • UN 1992 • FLAMMABLE LIQUID, TOXIC, N.O.S.(Aquagent® Titrant 2, titrant-component for Karl Fischer titration)
- IATA/ICAO: 3 II • UN 1992 • FLAMMABLE LIQUID, TOXIC, N.O.S.(Aquagent® Titrant 2, titrant-component for Karl Fischer titration)
- PAX: 305
- CAO: 307
- Store between 15°C and 25°C